

Angeltalk



www.angeltalk.com.cn


Tsinga Group, Flat Rm 1005 Wing

Yue Bldg 60-64 Des Voeux Rd West HK

TSINGA

Operation Manual

- ◆ Install the battery correctly, turn on the power knob, make sure the light on.
- ◆ Before listening, make the volume to medium to protect your ears from unpleasant voice.
- ◆ Put the earphones on, apply the gel and slowly move Angeltalk until you hear your baby's heartbeat.
- ◆ You can adjust the volume knob to catch the optimum signal.
- ◆ Fetal heart sounds is different from yours, fetal heart rate is 120-160bpm, its sounds like drumming. Please pay your patience to listen, you will get it.
- ◆ It's possible to hear sounds from 16 weeks into pregnancy.



Angeltalk

Notes:

Angeltalk prenatal heart listener brings you a clear sounds of your baby's heartbeat, but in order to protect your ears from unpleasant voice, we suggest you power on the device first before listening and make the volume in middle place.

Your baby and your body will make different sounds, at the time to listen, you can adjust the volume for the best effect.

You will be suggested to pay some patience to find the best position of fetal heartbeat.

Heartbeat sounds level has relations with fetal heart position. Sometimes, you may not hear because of the baby's movement.

In order to obtain a good sound effect, you're suggested to use Gel.

If not Gel, you also can put food oil or water.

Angeltalk probe should be touched your skin. Mobile phone, radio signals, etc will effect the Angeltalk use. Please keep them away you.

We suggest you not to detect fetal heartbeat more than 3 times a day, each time less than 3 minutes.

How to record the sounds

Recording to your PC:

- 1) Plug the recording cable into your PC and your Angeltalk Device
 - 2) Select : Programs
 - 3) Select : Accessories
 - 4) Select : Entertainment
 - 5) Select : Sound Recorder, and press Record
- Besides, you also can use other recording software to do it and make different style and format.

The best period of time to hear Fetal Heart Rate movement are:

1. Within 30 minutes after the pregnant woman gets up in the morning.
2. Within 60 minutes after the pregnant woman has lunch.
3. 30 minutes before pregnant woman goes to bed at night.

Trouble Shooting

If Angeltalk do not work, please check follows:

1. Volume knob whether or not put the lowest

volume.

2. Bad or inadequate battery power
3. Earphone cable whether or not good contact with the Angeltalk's Earphone jack.
4. Whether or not put the Gel
5. If you see the light on but the angeltalk was not working, that's lower battery caused please change a new battery.

Notes: it's a normal electricity sounds when you turn on the Angeltalk, not a bad or headphone's problems.

Disclaimer

Angeltalk is a non diagnostic device and does not replace the role of the doctor. You ought to consult the doctor about any pregnancy question. However, we still try our best to offer you the fetal heart rate for your reference. Do not as according as clinic diagnose but reference only during gestation. At the same time please operation it as instruction manual. Our industry and dealer will carefree if any of inadequately use.