


客户端下载体验 How to download App?



扫描二维码进入下载页面，或直接在 iPhone app store 搜索: E-touch。
Scan QR code below or search keyword: *E-touch* on the Apple App Store, you can download App.

下载安装手机App Download and install App on your iPhone or iPad.



将E-touch主机插入iPhone或iPad 音频口，将音量调至最大。将按摩贴片贴在需要按摩的部位，打开App。

Connect E-touch with iPhone or iPad through audio port, and maximize the volume. Stick two pasters on the body part and start the App.



选择模式，开始按摩 Select relaxing mode and start it.



名称(Name): E-touch
频率(Frequency): 9~45Hz
电流(Electric Current): <20 μ A
瞬间电压(Voltage): 130~350V
尺寸(Size): 58*26*14 (mm)
重量(Weight): 45g

温馨提示Warm Tips:

- 1、使用时，请从最低强度和最小频率开始逐渐提高；
During the operation, please gradually increase the strength and frequency from minimum.
- 2、请用湿毛巾擦净需按摩部位的皮肤表面，防止油污、化妆品以及其他污迹损坏按摩贴片的粘性；
Please keep your skin clean before sticking the pasters on it to avoid decreasing the pasters' stickiness.
- 3、不要在电源开着的时候用力挤压按摩贴，否则会受到强烈刺激；确定关闭电源后才能将按摩贴片移动到身体其他部位；
Do not squeeze the pasters when it is power-on. Ensure the device is power-off, before removing the pasters to other part of your body.
- 4、开车时禁止使用本产品；
It is prohibited during driving.
- 5、使用此产品过程中，如果感到不适，请立即停止使用；
Once you feel uncomfortable, please stop using it immediately.
- 6、在任何情况下，请勿拆卸主机；
Do not disassemble the device in any case.
- 7、使用完后，请先关闭电源，轻轻撕下贴片，并妥善保管；
After using, please turn off the power. Take the pasters off and keep the device well.
- 8、14岁以下儿童需在家长陪同下使用本产品；
Children under 14 years old should be accompanied by parents to use this device.



禁止使用者

It is prohibited for those people

心脏病或其他重大疾病患者
Patients with heart disease or other severe disease

使用心脏起搏器者
People fitted with a pacemaker

血压不正常者
People with abnormal blood pressure

孕妇或者产妇
Pregnant woman or puerpera

因淋浴或出汗后身体潮湿者
People with wet body after shower or sweating

关注我们

Please follow our WeChat:



E-touch 按摩体验从这里开始
Relaxing moment starts from here.

