

SMD Mini LED Phototherapy
User's Manual

Product & Accessory
Mini LED Phototherapy:

Like a "I-phone", it can be held by
hand, use with special gel
together.

Mini LED Device
Battery Charging(Can be
charged by computer with UBS
Connector)

User's Manual

Gift Packing

Recommended Special Gel:

For acne or problematic skin
care:

Acne Treatment Gel

For wrinkle& anti-aging skin
care:

Anti-aging& Skin Whiten
treatment Gel. Wrinkles Gel.

Mini LED Device's Function:

The Green Color LED reacts the
photosensitive effective
ingredient of Acne Treatment gel,
which reduces acne fungoid and
skin fungoid, reduces Oil
Secretion.

The Red Color LED reacts the
Photosensitive effective
ingredient of Anti-aging
Treatment Gel, Winkle Gel,
which make sure to improve
Regeneration of Skin, increase
Skin Elasticity, reduce Wrinkles,
make your skin bright and
elastic.

Usage:

① Battery Charging: 3 Hours is
enough, it mean the full charge
when indicating light is off.

Full charge can be used 75
Minutes, 5 times
(15mintues/treatment)

※ Please make the battery
charging before first usage.

② Clean the skin with clean
cream, then adjust the skin with
Toner, then open the pore with
warm water foam or hot-gas to
make sure the better absorbing
of Special Beauty Gel.

Spread the special beauty gel
on the face or necessary part by
massage method.

● Recommended Special Beauty
Gel

For acne or problematic skin
care:

AC treatment Gel

For wrinkle& anti-aging skin care:
Anti-aging treatment Wrinkle Gel

④ Keep 3-5cm between LED

device and Treatment Part,
make the treatment(Light
On) for treatment part till
the device is off
automatically (About
15Minutes is off
Automatically)

● For different skin care, choose
different LED Color.

Acne or Problematic skin: Choose
Green LED Color

Wrinkle/Anti-aging: Choose Red
LED Color

⑤ After treatment, please clean
the Special beauty Gel with
water or Cleaning Cream, then
use skin cream for repair and
adjust skin.

Treatment Session:

For Acne or Problematic Skin:

1 treatment/2-3days

For Wrinkle/Anti-aging skin:

1 treatment/day

Advantages:

1. For Acne or Problematic skin,
no pain, no side-effect. Low Cost

1 week can see the result.

Function:

- Removal Acne Fungus from
Dermis, restrain sebaceous
gland
- Removal keratin, Shrink
pores, let out sebum
smoothly.

For wrinkle/Anti-aging skin:

No Side-effect, low cost

2 weeks can see the result

Functions:

- Skin Rejuvenation and mitigate
wrinkles.
- Adjust skin situation and,
Shrink Pores