

Head & Eyes Massager

User Manual




Product Introduction

Head&Eyes Massager has been developed basing on principles of traditional Chinese Medicine Meridian treatment and with combination of modern science&technology. It starts with the Meridian in the head,stimulates the peripheral nerves,promotes the head microcirculation,opens channels and activates blood via massaging the head and effectively boosts the functions of the central nervous system to prevent cerebral arteriosclerosis,neurasthenia,headache and insomnia,stroke,dementia and sensory dysfunction of the central nervous system,movement disorders and other diseases,enhance memory,anti-aging,head to meet the day-to-day care of people and health care needs.


Parts name and function

Main body of Head Massager part



- 1、 Upper and lower adjustment knob
- 2、 Acupuncture point massage cone
- 3、 Airpipe
- 4、 Tightness Adjustment Button
- 5、 Windows for eyes

Controller



- 1. On/Off/Auto/Manual
- 2. "Heat" Button
- 3. "Mode" button
- 4. "Music" button
- 5. "Time" button for timing
- 6. Music mode display
- 7. Timing display
- 8. message program automatic display
- 9. message program manual display
- 10. Heating display

Product efficiencies

- ◎ Multi-directional Airbag Kneading&pressing acupuncture points in the head,help open the head Meridian,enhance blood circulation,remove eye fatigue and improve eyesight.
- ◎ Multi-point acupuncture Through the distribution of hundreds of biomass root cones to acupuncture on the head to increase the oxygen exchange function of the brain,eliminate head pain and improve head health.
- ◎ Comfortable brain wave music to make you sleepy With the effect of the natural music, remove the tension of the cerebral nerves, release the spiritual stress, improve the self-recovery of the brain cells, stimulate the potential of thinking.
- ◎ Gentle Heating Function Soft warm effect can not only promote the blood circulation, but also reduce the pain and sense of tension in head & eyes.


User Guide

◎ Preparations before use

- 1、 Open the battery cover on the back of the controller.
- 2、 Load four AA(5#) batteries correctly according to their positive and negative poles.
- 3、 Place back the cover after the battery loading.

◎ Put on and Adjust

- 1、 Put on the massager around your head,adjust the knob to make sure the eyes can be fully covered.
- 2、 Use the buttons at both sides to adjust to the comfortable tightness,according to the displayed digits at both sides to get the balance,a little tight will be better.



◎ Use and Operate

【Auto Mode】

- 1、 Press the button numbered 1 above(on/off) to enter Automatic status,the related light to indicate Auto will be on automatically.
- 2、 Under the Auto Status, buttons numbered 2 (||||),3(✕),4(♪) and 5(⊖) will have changes accordingly with different light indicator,please refer to Operation Directions for every single function details.

【Manual Mode】

- 1、 Press the button 1 above,you can enter Auto mode,press again to Manual mode.(will appear by turns between the two modes by pressing this button)
- 2、 After entering the Manual status,only the music can be heard,and the others stand-by,you please operate on button 2,3,4,5 to get the corresponding functions,the light indicators will change to display accordingly.

【"Mode" button】

Five intelligent massage programs are available and switchable by pressing on Button numbered 3 above.

Mode for rest and peace

Via airbag to knead and press alternatively,under the guidance of the natural music,head fatigue and tension can be removed and you will be led back to nature and easy status of moods.

Mode to activate your thinking.

Through the strength massage to Baihui,Taiyang acupoints,as well as auxiliary massage and vibration to other acupoints,together with warm temperature to comfort,boost head blood circulation,enhance vitality of brain,evidently improve thinking ability and memory of your brain!

Mode for Natural Sleeping

Via comfortable music of hypnosis and gentle,smooth airbag massage,guided the brain to dormant state,effectively release tension,improve sleeping quality.

Mode to activate the brain and comfort your heart

Via the gradually stronger and stronger kneading,press and multi-frequency vibration to improve the head microcirculation,enhance the central nervous system,prevent and treat cerebral arteriosclerosis,stroke and other brain diseases.

【Music Button】

Eight natural light wave music are available for this massage and switchable by pressing on Button numbered 4 above.

● Sound of Basi Flowing water

The flowing water is just like the surging of life, calmly washes away the noise and glitz in the world,the listeners can obtain direct inner harmony and tranquility from it.

● Basi Sea wave sound

With the fluctuations and soothing of the music,the listeners will feel staying over the sea,the left and right brain can communicate more smoothly,the door of wisdom can be gradually opened accordingly.

● Basi Chirp Sound&Music

The perfect ensemble of the cheerful sound of birds and peaceful music has formed the greatest masterpiece of the nature,and the listeners can enter the spiritual journey of Heaven & Man into one body.

● Sound of the Peak Rain

The dry soul is just like the cracked earth and needs the moisture of the rain.To listen to the sound of rain at the depth of nature can make your heart satisfied and spirit harmonious.

● Peak Sound of night,Sound of all animals

You can feel Secluded Orchids from the vacant valley,Chirps from the deep gullies,together with the response sound from other animals,the lives of the Big Nature can bring mental & physical happiness to human beings and we can see vividness from everything in the world.

● Blue Bubble sound

A breath from the depth of the earth,comes out from the water surface,the natural breaths can adjust the hormone within the body,enhance metabolism and improve bodily functions.

● Sound of the water at night in the forest

Pure sound of night,flowing water can gently wash away the chaos in the mind from the listeners,remove tension and release the psychological stress.

● Sound of birds in the forest

Attracted by the melodious and cheerful chirping sound of the bird,the listeners can temporarily go far away from the uproar of noisy society,fresh smartness can be infused into the tired and sleepy brain.

【Timing Button】

Five different time periods are available for timing:5,10,15minutes respectively.The default value is 15m,you can set as you like and the time chosen will display on the remote control.

【"Heat" Button】

Three different temperatures can be provided,the indicator will change accordingly with each different press.

- 1、 Low(The indicator will glitter once.)
- 2、 Middle(The indicator will glitter once.)
- 3、 High(The indicator will glitter once.)

【Switch off】

Press the number 1 button(on/off) for 2 secs or so to switch off the operation.



Note: Please wait 10 seconds to reset for the next use!

Specifications

Product Name: Head & Eyes Profession Massagist
 Product Model:
 Operational Environment: Minus 30 to plus 55(Celsius degree)(Temperature)
 Vibration Frequency: 40-100Hz
 Air Pressure: 86Kpa-106Kpa
 Operational Voltage: 6V(4xAlkaline AA#5 battery)
 DC Adapter: Input 220-240V 50/60Hz 0.5A
 Output 6V== 1500mA
 Power Consumption: Less than 6W
 Dimensions: 145(L)×68.5 (W) ×45 (H) mm(Remote Control)
 230(L)×220 (W) ×234 (H) mm(Massager)
 Net weight: 0.72 KG
 Attestation: CE, RoHs

Attentions

▲ Persons listed below please don't use this massager.

- 1、 With head or eyes disease or once wounded but not recovered yet.
- 2、 With swelling near the head caused by heavy knock or hit.
- 3、 With poor health or serious patients still under treatment.
- 4、 Woman in pregnancy or physiological period.
- 5、 Persons who are using internal electronic medical equipments such as Pacemaker, etc.
- 6、 Children under 12.

▲ Safety Tips for use

- 1、 Not suitable to be placed or used under high temperature or humidity surroundings.
- 2、 Please don't put towel or else to cover the massage when using to avoid excess heat accumulation.
- 3、 Please quit using the massager immediately if abnormal feelings or severe pains occur.
- 4、 Please control the using time to be within 15minutes each time.
- 5、 Please don't use after drinking or when you don't feel well.
- 6、 Persons with mental incompetence or physical dysfunctions, please use under the custody of others.
- 7、 To make best use of this massager, please choose Alkaline batteries and maintain adequate power.
- 8、 Please don't lie down or walk or drive outside when using this equipment.
- 9、 Please try not to drop the apparatus onto the floor or immerse into the water.
- 10、 Please don't dismantle the apparatus of its plastic frames, so as not to lead to electric shock or other parts to malfunction.

Service Guarantee

For this product since the date of purchase, we guarantee to provide the following services:

1. For failures under normal operations, we offer free repair and free exchange of defective parts.
2. We provide service only but not for free for the following conditions:
 - a. Warranty expired.
 - b. Failure caused by wrong operation or unauthorized demolition.
 - c. unable to show this Guarantee.
3. Maintenance Details:
 - a. If faults occur under normal operation according to instructions during warranty period, please send the product to the local dealer or send to our company for maintenance directly.
 - b. For failures caused by human, we can offer free maintenance, while all accessories and freight at user's cost.
 - c. For maintenance outside the warranty period, free maintenance will be available, while all accessories and freight at user's cost.

Warranty Card

Product Name:	Head & Eyes Profession Massagist
Product Model:	
Tester:	
Production date:	
Name of Buyer:	
Buyer's address:	
Buyer's telephone:	
Zip code:	
Purchasing date:	
Purchasing place:	