

Operating instructions of the multifunctional electromagnetic belt

MRK-505 operating instructions

- ※ The machine is designed with six buttons [upper left is intensity button "-", the middle is ON/OFF buttons(stratup and shutdown), upper right is intensity button "+". Lower left is HEAT button(sauna), the middle is EMS(low- frequency pulse), lower right is VIB button(vibration)]
- ※ On starting up, LCD shows the whole state, what display screen shows is the default mode of strating up: MODE : 1、 LEVEL : 00 ; VIB : 0 ; HEAT : 0
- ※ The view of the main part



- ※ Accessories: an operating manual, a volume about health knowledge, a power adapter dedicated for the machine.
- ※ Time : working time of the program default is 15 minutes
- ※ Operation instructions of buttons
 1. Connect the dedicated power adapter, press the button Φ , start up after hearing the sound " beep ". After startup, enter automatically the default state : the backlight of the LCD screen is on, the red LED lights, the intensity of the heating function(HEAT) is "0", discharge function(EMS) is for mode (MODE) intensity is "00", intensity of vibration function (VIB) is "0".
 2. After startup, in any time without any operation or function output, in this case, within 20 seconds the machine shut down automatically.
 3. Heating (HEAT) function: each time you press the HEAT button, the intensity of heating will circulate in this way "0—1—2—3—0"
 4. Low frequency discharge(EMS) function: each time you press the button EMS, low-frenquency discharge will convert a pattern, by "1—2—3—4" these four model cycles. Press the button "+", you can increase the intensity

of low-frequency discharge, each time you press, the intensity increases "1", and "10" is maximum, at this point to the maximum, and then pressing the button "+" is invalid. Press the button "-", you can reduce the intensity of the low-frequency discharge, each time you press, intensity reduces "1", and "00" is minimum, so far as the minimum, then pressing the button "-" is invalid. (Note: every conversion model, the intensity of the new pattern begins from "00".)

5. Vibration (VIB) function: each time you press the button VIB, vibration intensity will circulate from "0-1-2-3-0".
 — Note 1: 3), 4), 5) these three functions, can be executed simultaneously, also can be used alone.
 — Note 2: when using the "low-frequency discharge" funcnion, bipolar electrode slice must cling to the skin to work properly, in case the sudden current enhance wounding, we made such a technological protection measures.
6. (This is new). Increase the temperature control function to control the temperature within a certain range (such as: surface temperature will be locked up to 65 degrees). This function is especially important in high heat grade (HEAT-3)

※ Technical parameters

- 1, the normal working conditions
 - a) Ambient temperature: 10 °C -40 °C
 - b) Relative humidity: <80%
 - c) Atmospheric pressure: 700hPa-1060hPa
 - d) Power adapter :110-220V 50Hz, 5V / 1.5A dedicated power supply
- 2、 Performance index
 - e) Timing range: 15 min
 - f) Power consumption: <0.5W
 - g) Skin electrode low-frequency pulse output frequency / S: 4 times, 8 times, 15 times, 30 times
 - h) Output waveform: square wave, rectangular wave, sawtooth, triangle wave, peaked-wave, step wave.
 - l) The heat treatments temperature: 30 degrees, 40 degrees, 50 degrees, 60 degrees, 70 degrees.

※ Common faults and solutions

Fault phenomenon	Cause of the fault	Solution
Can't boot up or crash	The power isn't connected well	plug in the charger and restart
No feeling of pulse	check whether the belt is close to the skin, and if the intensity of pulse has been increased	Fasten tightly the belt with the skin, and adjust the pulse by pressing the buttons "+", "-"
No feeling of sauna	check whether the belt is close to the skin, and if the intensity of sauna has been increased	Fasten tightly the belt with the skin, and adjust the HEAT(sauna) button
No feeling of fat rejection	check whether the belt is close to the skin, and if the intensity of fat rejection has been increased	Fasten tightly the belt with the skin, and adjust the VIB(vibration) button

※ After-sales service for the series of Merkon Brand, warm tips as follow:

★ Counseling service on telephone

Announcements

To ensure this kind of massage brings you happiness, the following is prohibited to do. Or used under the guidance of professional doctors.

1. Please do not use it when driving.
2. Please do not use it beside the interferential equipment with high temperature, or damp places.
3. Please use the dedicated DC 5.0V power adapter configured by our company, and pay attention to electrical safety.
4. Body has congenital diseases, red and swollen, inflammation, congestion, consumptive disease, and other symptoms, those who have had abdominal surgery and have not been healed, for those persons, should under the guidance of professional doctors when using it.
5. People who suffer from high blood pressure, heart disease or other diseases, should under the guidance of doctors using this product.
6. It's forbidden to use it for those people like under 20 years old, mental patients, the old and infirm.
7. This product is just for personal health care, and for no other purpose.

After-sales service and maintenance warranty service

Maintenance

1. Take measures to anti-moisture, anti-rats, anti-pollution. Keep away from high heat, sources of ignition zone; avoid glare.
2. Clean the handle controller should be performed in the state of shutdown. When cleaning the product, you should wring the wet towel to wipe. In case of stubborn stains, can be dipped in a little alcohol wipe.
3. After the use of the product, it should be placed in a safe place, do not hold heavy objects on the product.
4. When a fault occurs, maintenance of professional and technical person should be invited, do not disassemble products.

3

Commitment

From the date of sale, for performance failure, in seven days, with product warranty card, certificate and invoice, you can go to the original purchaser counters handle retirement / replacement procedures, return / exchange of the products must ensure that there is no artificially external damage, power adapter, packaging, manuals and other accessories are intact; Since the date of sale, in 15 days, if performance failure appears, you can choose to change another product in the same price or more expensive or you can ask for repair. Please contact sales department.

The following does not belong to free warranty, but can be provided with maintenance services.

(Free maintenance, charge fees only for accessories) for accessories replacement you should ask for the price before the payment

1. Due to customer's unsuitable use, storage, maintenance, damages appear.
2. Damage caused by repairing and dismantling without authorization
3. Without our products warranty card and purchase invoices.
4. Damage caused by force majeure.

Warranty methods

This product can enjoy one year free warranty from the date of purchase, and can enjoy nationwide assurance of maintenance. When repairing this product, please show repair warranty card, purchase invoice (free repair promised by the company includes only production time and material costs, transportation costs are not included, and the company does not provide on-site maintenance services). We are willing to assume relevant obligations such as recycling, treatment, re-use of the waste products.

Tips: Final explanation of this product belongs to Merkon, for the upgrading of the product, excuse us that we will not give notice in addition!

Warranty card:

With this warranty card can be maintained Date:

(without this card can't be maintained)

User name	Sex	Age
Address		
Telephone	Postcode	
Purchase place		
Purchase price	Sales order number	
Trade name and specification	Serial number	
Dealer signature	User signature	

4