


Beauty Care Machine

Operation manual

**Thank you very much
for your purchase**

- Thank you for choosing I-PH-S beauty care machine
- Please read this operation instructions manual carefully before use and keep this brochure for future reference
- This product is designed for beautifying skin, don't use or sell it as medical purpose.



Contents

Specification	1
Safety Precautions	3
Important notice	4
Quick user guide	9
Warranty card	12

Specification

Item No.: FF0382

Name: Ionic-Photon-Ultrasonic Skin Care Machine

Adaptor input power: AC 100-240V

Adaptor output power: 15V, 400mA

Dimension: 17.8×4.8×4.5cm






Weight: 175g

Ultrasonic frequency: 3MHz±5%

Foreword

Skin is the largest organ of human body, related to the health and image. More and more people, especially women, take more attention to the skin. However, gravity, sunshine, radiation, etc. are harmful to skin, causing sag, loose, wrinkles, spots, acne and other skin problems. To protect and beautify skin, we need not only skin care products, but also a professional beauty device to get better effect.

Home use I-PH-S beauty care machine combines three popular & effective skin care technology: Ion (galvanic current), photon (red, blue, green lights), ultrasound. Use with skin care cosmetics, it's convenient to care skin at home, easy to carry and operate. The efficacy is as below:

-  Clean skin deeply
-  Eliminate wrinkles, smooth skin
-  Whiten & tighten skin
-  Reduce fat, shape slim face contour
-  Dilute spots, dark circles, eye bags

Safety Precaution

Note: Please read following information carefully before using your machine. This will help you to ensure the warranty validity and extend the life cycle of this product.

1. Store the item away from the humidity and do not use it under water, otherwise the electronic components maybe corroded or get short circuit.
2. Never attempt to disassemble the product by yourself. All the components can be only maintained and repaired by professional technicians.
3. Don't use any other except the original one. Using incompatible may cause malfunction or danger.
4. Any people with the following conditions must avoid using this product: patients with malignant; women during her pregnancy or period; patients with dermatosis of foci on hemorrhage; urgent patient with unidentified diseases.
5. Keep this product away from children.

Important Notice

1. This machine should be used under normal condition. Remember to turn off the power and unplug the cable when you stop treatment. This will avoid over-heating problem caused by the air vibration and to prevent any to shorten product life cycle.
2. It's normal that the massage head will be a little hot when using ultrasound, and will feel pricking. Using skin care product could relief the discomfort of heat. Gels are recommended.
3. Avoid using this machine in one area for too long time. Move it around on your skin to avoid focusing effect caused by the wave energy.
4. In order to get better treatment, attach the massage head to skin completely when using.
5. Please clean up the massage head after use. A wet soft cloth is recommended to prevent the scratch. Don't use any chemical cleanser to clean the product, otherwise it maybe causes malfunction or shorten product's life cycle.

The Simple Allergy Test:

1. Clean the wrist.
2. Apply some essence that you intend to use on the cleansed area, let it dry naturally (During the test, please don't touch water nor wash the essence away). If you experience allergy or any skin abnormalities, stop the test immediately and contact your doctor or nearest hospital straight away.
3. After 48 hours, if there is no allergy and abnormality, you may use this device to enjoy the beauty care process at home.



Nomenclature

Dynamic lights:

- ▼ 12pcs of red LED: 630nm
 - ▼ 12pcs of blue LED: 465nm
 - ▼ 12pcs of green LED: 520nm
- High energy, highly focus

High quality antioxidant stainless steel massage head



High-grade galvanized sidebar



LCD Display

Sonic)))))
Ion + @@@@
Ion - @@@@
Photon ☀
Time 15 Min
Bring Beauty Home

Control Panel



Power light function:

1. Red light: 630nm

It can reach the dermis layer to stimulate fiber cells, improve blood circulation, activate skin, reduce wrinkles and smooth skin. Suitable for any skin type.

2. Blue light: 465nm

It can kill bacteria, calm and against skin rashes & acne, strengthen the protein fiber organization, shrink pores. Suitable for oily skin type.

3. Green light: 520nm

It can relax skin and relieve the pain of skin. Suitable for the mixed skin type.

Operation Notice


- © Plug the adaptor into device and the socket, there will be a sound "Beep" and the machine is ready to work.
- © Apply gel or essence on face or massage head, switch on the machine.
- © Attach the surface of massage head on the skin completely to move according to the indicated graphics.
- © If feel too dry to move, add some more gel or essence to continue on. (Dosage: 4-5 g per time)

Basic way of massage:

Lift bottom-up

Left and right symmetrically

Slow movement



OPERATION

Function combination	Sonic	ION+	ION+	Photon
Single function	Orange arrow up	Blue arrow up	Red arrow up	Green arrow up
Double functions	Orange, Blue, Red, Green arrows up	Blue, Red, Green arrows up	Red, Green arrows up	Green, Blue, Red arrows up
Three-function combination	Orange, Red, Green arrows up	Orange, Red arrows up	Red, Green arrows up	Orange, Red arrows up

1. Switch on the machine → LCD light up after "Beep" sound.

2. Single function:

Press "MODE" to select function, say, "Sonic" → "Sonic" flashes → press "SHIFT" to adjust power → Press "CFM" to work.

3. Mixed functions (2 kinds of function):

After single function proceeding → Press "MODE" to choose the other function you need, say, "Ion-" → Ion- flashes → Press "Shift" to adjust power → Press "CFM" to work.

4. Mixed functions (3 kind of functions):

After 2 functions proceeding → Press "MODE" to choose another function you need, say, "Photon" → "Photon" flashes → Press "SHIFT" to choose LEDS → Press "CFM" to work.

5. Timing

The default setting of auto off is 15 minutes → Press "MODE" to choose "Time" → "Time" flashes → Press "SHIFT" to adjust Time from 1-19 mins Press "CFM" to work.

Remark: Ion+ & Ion- can not work at the same time.

Tip: Simple Testing for Ultrasound

First step: Drop a little water on the top of massage head

Second step: Switch on the machine, select mode "Sonic" and power, press "CFM" to start working.

Third step: Water splashing means ultrasound works well.

Restart the machine to continue using after auto shutdown.



Quick User Guide

Daily care. (Please don't use the device more than 3 times within one week and don't exceed 30 minutes each time.)

First step: Clean skin deeply

1. Thoroughly remove make-up and impurities first.
2. Apply proper cleanser on the face.
3. Select ultrasonic plus ion+ function with comfortable intensity.
4. Along the indicated route, move the massage head from jaw to the bottom of ear slowly. Other place as the indicated graphic.
5. Clean face with water.

Second step: Lead in nutrient and shape face contour.

1. Apply proper gel or essence on the surface of massage head.
2. Select ultrasonic plus ion- function with comfortable intensity.
3. Along the indicated route, move the massage head from jaw to the bottom of ear slowly. Other place as the indicated graphic.

Third step: Photon care

1. Select a light as demand (red/blue/green)
2. Illuminate light on the skin where need care. Move the massage head every 5 seconds.

★ To shorten the treatment time, the second step and third step can work together.


Dilute spots & whitening (one time every two days)

First step: Clean skin deeply (please refer to the first step of daily care)



Second step: Lead-in whitening essence

1. Switch on the machine and select ion- function with comfortable intensity.
2. Apply proper whitening essence (such as VC) on the spots or skin where need whitening, use machine to work on it. Let the massage head to stay about 20-30 seconds in one position then move to the other.



Acne treatment (One time every night)

First step: Clean skin deeply (Please refer to the first step of daily care)

Second step: Lymphatic drainage

1. Switch on the machine and select sonic function with comfortable intensity. perform lymphatic drainage according to the indicated route.
2. Massage around acne.
3. Use blue light to illuminate on the acne area about 15-20 minutes
4. Use red light to care about 10 more minutes.



Clean the skin



Lymphatic Drainage

Sensitive skin (one time every 2 days)

1. Clean skin deeply (please refer to the first step of daily care)
2. Switch on the machine and select sonic function with comfortable intensity.
3. Use green light to illuminate on skin where needing care about 15-20 minutes.



Acne inflammation

Tip:

User should touch the galvanized sidebar when using ion+ & ion- functions.



[Guarantee Certificate]

Product model:_____ Purchase date:_____

Purchase Place:_____ Sales Dept:_____

One year Repair Guarantee

In the unlikely event that your product develops a fault within the first one year of purchase: Please return your product along with your receipt and the original guarantee certificate and we will repair it free of charge provided you have operated it according to the operating Instructions

What is not covered?

1. Breakdown or damages due to
 - (1) not following the operating instructions.
 - (2) incorrect or unauthorized repair
 - (3) accidental damages during your own transportation, or you drop it on the floor or you deliberately damage or neglect the product after purchase.
 - (4) fire, earthquake, weather such as flood, lightning, high wind.
2. Please keep your guarantee certificate safe, as no duplicated certificate will be issued.

2015007000:

Date of production:

Certificate

Model: FF0382

Checker: